

Arrow Epidural Products

CLINICIAN-INSPIRED PRODUCTS FOR REGIONAL ANESTHESIA



SJ-05501
Arrow FlexTip Plus Kit

ARROW FlexTip Plus® EPIDURAL CATHETERS

Designed to reduce complications in the epidural space. Excellent kink-resistance and tensile strength with a soft, flexible tip reduce the risk of paresthesia and inadvertent penetration of blood vessels or the dura, which can misdirect delivery of anesthetics.^{1,2} Simpler and faster epidural procedures are facilitated by the catheter's unique coil reinforced body. The structural integrity of this design provides resistance to collapsing and requires no stylet for insertion. Catheter radiopacity is enhanced by the inner spiral-wound stainless steel coil body construction.



ARROW TheraCath® EPIDURAL CATHETERS

Designed for pain management procedures where increased steerability is an important factor. The kink-resistant catheter provides high compression strength to resist collapsing. Smooth, durable polymer exterior and inner spiral-wound stainless steel helix provide column strength for directionability and enhanced radiopacity. Removable stylet further enhances positive placement. Unique omniport tip design provides uniform disbursement of anesthetic.



1 Junega M, MD, Kargas GA, MD, Miller D, CRNA, Perry E, MD, Gupta B, MD, Pajel V, MD, Botic Z, MD, and Rigor B, MD. Comparisons of Epidural Catheters Induced Parasthesia in Parturients. Poster presentation from the American Society of Regional Anesthesia (A.S.R.A.) meeting: presented April, 1995 – Orlando, FL. Regional Anesthesia 1995; 20 (Suppl.):4.

2 Junega M, MD, Kargas GA, MD, Miller D, CRNA, Perry E, MD, Botic Z, MD, and Rigor B, MD. Incidence of Epidural Vein Cannulation in Parturients with Three Different Epidural Catheters. Regional Anesthesia 1996; 21 (Suppl.):4.

For ordering information and a list of Arrow Epidural products see the back of this sheet.

REGIONAL ANESTHESIA

ARROW
INTERNATIONAL

Products with FlexTip Plus® (19 Ga.) Catheter

Product No.	Key Components
AK-05500	17 Ga. Sharp Tuohy; 5cc glass LL LOR syringe; 10/cs
AK-05501	17 Ga. Sharp Tuohy; 10cc plastic LS LOR syringe; Clear Drape; 10/cs
AK-05502	17 Ga. Sharp Tuohy; 5cc glass LS LOR syringe; Clear Drape; 10/cs
AK-05503	17 Ga. Dull Hustead; Attachable wings; 5cc glass LS LOR syringe; Clear Drape; 10/cs
AK-05503-L	17 Ga. Tuohy; Attachable wings; 5cc glass LL LOR syringe; Clear Drape; 10/cs
AM-05500	17 Ga. Sharp Tuohy; Advancer; 5cc glass LL LOR syringe; 10/cs
BP-05501	17 Ga. Sharp Tuohy; Advancer; 10cc glass LL LOR syringe; Clear Drape; 10 /cs
EC-05400-E	17 Ga. Dull Hustead; Attachable wings; 0.2 µ filter; Set; 20 /cs
EC-05500	FlexTip Plus Epidural Catheter with SnapLock only; 20/cs
JC-05400-E	17 Ga. Dull Hustead; Attachable wings; 10cc plastic LS LOR syringe; Set; 20/cs
SJ-05501	17 Ga. Sharp Tuohy; Advancer; 5cc glass LS LOR syringe; Clear Drape; Epidural Catheterization Kit with full prep tray; 10/cs
SJ-05501-0	17 Ga. Sharp Tuohy; Advancer; 5cc glass LS LOR syringe; Clear Drape; Epidural Catheterization Kit with basic prep tray; 10/cs
SM-05501	17 Ga. Dull Hustead; Attachable wings; Advancer; 5cc glass LS LOR syringe; 10/cs
WJ-05401	17 Ga. Dull Hustead; Attachable wings; Advancer; 5cc glass LS LOR syringe; 30" x 24" Clear Drape; 10/cs

Products with FlexTip Plus (20 Ga.) Catheter

Product No.	Key Components
EC-05420-P	18 Ga. Super Dull Crawford; Advancer; 10cc plastic LS LOR syringe; 0.2 µ flat filter; Caudal Placement Set; 20/cs
EC-05520-P	18 Ga. Dull Tuohy; Advancer; 10cc plastic LS LOR syringe; 0.2 µ flat filter; Lumbar Placement Set; 20/cs
AK-05420-P	18 Ga. Super Dull Crawford; Advancer; 18" X 30" Clear drape; 5cc glass LS LOR syringe; 0.2 µ flat filter; Opsite® Dressing; Caudal Placement Kit; 10/cs
AK-05520-P	18 Ga. Dull Tuohy; 18" X 30" Clear drape; 5cc glass LS LOR syringe; 0.2 µ flat filter; Opsite® Dressing; Lumbar Placement Kit; 10/cs

Products with TheraCath® (19 Ga.) Catheter

Product No.	Key Components
AK-05000	17 Ga. needle; Attachable Wings; 5cc glass LL LOR syringe; 10/cs
EC-05000	19 Ga. TheraCath Epidural Catheter set; with SnapLock; 20/cs
MM-05001	17 Ga. Dull, Hustead; Attachable wings; 5cc glass LS LOR syringe; 10/cs
SL-05000	17 Ga. Tuohy; 5cc glass LS LOR syringe; Clear Drape; 10/cs



AK-05000
Arrow TheraCath Kit

Combined Spinal Anesthesia / Epidural Catheter Products

Product No.	Key Components
AK-05560	17 Ga. Sharp, Tuohy; 26 Ga. Gertie Marx® needle; 19 Ga. FlexTip Plus catheter; 5cc glass LS LOR syringe; Clear Drape; 10/cs
UI-05502	17 Ga. Sharp, Tuohy; 27 Ga. Whitacre® Pencil Point Tip needle; 19 Ga. FlexTip Plus catheter; 5cc glass LS LOR syringe; Clear Drape; 10/cs

Single Shot Epidural Products

Product No.	Key Components
AK-03000	20 Ga. needle; 5 cc glass LS LOR syringe; 17 Ga. X 1" LS introducer; 20/cs
SA-03001	18 Ga. Dull Tuohy; 5cc glass LS LOR Syringe; Clear Drape; 20/cs
SH-03001	17 Ga. Dull Hustead; Attachable wings; 5cc glass LS LOR Syringe; Clear Drape; 20/cs
SL-03000-1	Epidural Access Tray; Contains procedural injection needles, syringes, drape and prep supplies; 10/cs

Epidural Accessories

Product No.	Key Components
AN-02017	17 Ga. Continuous Epidural Needle; Dull Hustead; 100/cs
AN-05501	17 Ga. Continuous Epidural Needle; Sharp Tuohy; 100/cs
AN-05505	17 Ga. Continuous Epidural Needle; Tuohy; 25/cs
AS-05500	SnapLock Catheter / Syringe Adapter; 50/cs
EF-05500	.20 µ Flat Epidural Filter; 100/cs
LR-05501	10 cc Plastic LS LOR Syringe; 25/cs
LR-05502	5 cc Glass LS LOR Syringe; 25/cs
LR-05503	5 cc Glass LL LOR Syringe; 25/cs

Caution: U.S. Federal Law limits this device to sale by or on order of a physician. Contents of unopened, undamaged package are sterile. Disposable. Refer to package insert for current warnings, indications, contraindications, precautions, and instructions for use.

All Arrow Epidural products and accessories are Latex-Free.

Gertie Marx is a registered trademark of Walter A. Zohmann.

This is a listing of the most popular Arrow Epidural products. For additional kit configurations, contact your Arrow Sales Representative or Arrow Customer Service.

To locate an Arrow Distributor or for ordering information, contact Arrow Customer Service at:

United States

Tel: 800.523.8446

Fax: 800.343.2935

www.arrowintl.com

International

Tel: +1.610.655.8522

Fax: +1.610.655.8566

Distribution Worldwide:

Arrow offices are located in Belgium, Canada, Czech Republic, France, Germany, Greece, India, Italy, Japan, Mexico, Netherlands, Slovakia, South Africa, Spain, United Kingdom and the United States.

Arrow International, Inc., 2400 Bernville Road,
Reading, PA 19605 U.S.A. Tel: 610.378.0131

www.arrowintl.com

EPD-S 04/2007 CRB-07-0101

© 2007 Arrow International, Inc. U.S.A. All rights reserved.
Printed in the U.S.A.

